



CLEANSIFY

SMART WIRELESS EAR CLEANING OTOSCOPE



User Manual

Package Contents :

- 1 Metal and silicone ear spoon x 1
- 2 Q-Elastic Ear Spoons x 1
- 3 Charging Cable x 1

Download the APP


- IOS
- Search for "bebird" in App Store
- Android
- Search for "bebird" in Google Play



OR scan the QR code to download the app





- 2 Open bebird APP select bebird R1
- 3 Press and hold the power button to turn on the ear cleaning rod




Connect Product WiFi

- 1 Connect product WiFi "bebird-R1-XXXXXX"
- 2 Back to bebird APP Click on "Start" to use it



Replace the ear spoon head

To remove ear tip twist anticlockwise.
To fix the replaceable ear tip twist clockwise



Insert and fasten the ear spoon head

The fixed ring must be tightened to stop the picker's falling.

To avoid damage caused by material aging, please do not use the same ear tip for a long time. We advise you to replace the ear tip at least once every three months

Product parameters

Pixel	300W
Lens diameter	3.5mm
Optimum focal length	12 ~ 15mm
Depth of field	10 ~ 50mm
Image transmission rate	20fps
Working temperature	-10 ~ 45°C
Gyroscope	3axis
Battery	180mAh
Battery life	1h
Charging time	1h
Input current	5V = 0.5A
Rod weight	13g
Working frequency	2.4Ghz
Network standard	IEEE 802.11b/g

Troubleshooting

Problem	Solution
No response after startup	Please use the incidental charging cable to charge. The charging time is one hour
The camera doesn't turn on	Cover the cap for a few seconds and remove it to restart. Try a few more times
The app exits quickly	Delete the existing app and scan QR to reinstall it
Feel tricky with replacing ear tips	Align the replacement ear tip to camera and tighten clockwise
Indicator light flashes	Please use the incidental charging cable to charge. The charging time is one hour
Abnormal display in app	Interrupt the connection, restart the product to connect again
Unstable Wi-Fi connection, disconnected occasionally	Cover the tap for a few seconds and remove it to restart, reopen the app to connect again

If you still have problems with the above procedures, please contact our Customer Service Representatives directly for further assistance

Attention

- 1 When cleaning lens, please wipe it carefully with professional alcohol cotton swab.
- 2 Please don't use it while you are moving, and protect from being hit by others.
- 3 Please do not put the product in liquid to avoid damage.
- 4 Product built-in rechargeable lithium battery, so it must be charged once a month If it is not used frequently.
- 5 Avoid product exposure to sunshine directly, especially the accessories, in order to avoid softening.
- 6 During using it, the temperature of the product will rise slightly (up to 35 degree), please rest assured to use.

www.portronics.com

For Support

Visit : www.portronics.com

Email : help@portronics.com

Contact No. : +91 9555245245 , 01-42413131

(Working Hrs. : Monday - Saturday 10.00am - 6.00pm)



For E-waste & ROHS compliance refer to <http://www.portronics.com/ewaste.html>

