Instruction Manual

Electronic Blood Pressure Monitor

Model: BP-205

- EN** User manual in English - 1
- DE** Anleitung auf Deutsch - 16
- IT** Istruzioni in italiano - 32
- FR** Instructions en Français - 48
- ES** Instrucciones en Español - 64

File No.: BRAV_BP_011
Version: 1.3
Revise date: 2023-07
Made in China

CONTENTS

- Safety Instructions 2
- Working Principle 4
- Intended Use 4
- Product Introduction 4
- Product Structures 4
- Symbol Explanation 5
- Product Specification 6
- Settings 8
- Blood Pressure Measurement Method 8
- Error and Fault Troubleshooting 11
- Maintenance and Cleaning 12
- Battery Replacement 12
- After-sale Service 13
- EMC Information and Statement 13
- Manufacturer Information 15

1

1. Safety Instructions

Thank you for purchasing Runstar Arm-type Electronic Blood Pressure Monitor.
⚠ Please read this instruction manual carefully before use.
⚠ This device is not a direct diagnostic device. It is only suitable for monitoring the blood pressure and pulse rate of adults. The measurement results should not be used for self-diagnosis or treatment. Please follow your doctor's advice.

Blood pressure is an important indicator of the human circulatory system. The home blood pressure monitor is easy to use regularly. It is crucial for preventing cardiovascular disease, encouraging hypertension patients to seek treatment, and managing high blood pressure.

Blood pressure fluctuates continuously throughout the day and night. The highest values typically occur during the day and the lowest at midnight. It is recommended that you measure your blood pressure at roughly the same time each day.

Various activities such as physical exercise, excitement, stress, eating, drinking, smoking, body posture, and other factors (including taking a blood pressure measurement) can influence blood pressure readings. Due to these factors, it is uncommon to obtain identical blood pressure readings in multiple measurements.

Definition and Classification of Blood Pressure Levels from WHO-ISH

Classification	Systolic Pressure (mmHg)	Diastolic Pressure (mmHg)
Ideal Blood Pressure	<120	<80
Normal Blood Pressure	<130	<85
High Normal	130~139	85~89
Hypertension Grade I (Slight)	140~159	90~99
Hypertension Grade II (Medium)	160~179	100~109
Hypertension Grade III (Serious)	≥ 180	≥ 110
Isolated Systolic Hypertension	≥ 140	<90

The information above is for reference only.

1.1 Warning:

This product is intended for daily blood pressure monitoring only and should not be used for diagnosing hypertension.

- Use this product only for blood pressure measurements. Using it for other purposes may lead to accidents.
- This device is suitable for use by adults only. Do not use it on infants, toddlers, children or individuals who cannot express themselves.

2

- Self-diagnosis and treatment based on blood pressure readings are not recommended. Consult healthcare professionals for proper evaluation and guidance.
- Patients with severe blood pressure or circulation disorders, as well as blood diseases, should use this device under medical supervision.
- The displayed pulse reading should not be used as medical heart rate monitoring.
- Patients with arrhythmia or atherosclerosis should have their blood pressure readings confirmed by medical professionals.
- Do not apply the cuff on limbs with intravascular access, therapy, or arteriovenous shunts to prevent blood flow interference and patient injury.
- Do not use the cuff on an injured or operated arm to prevent further harm.
- In the presence of common arrhythmias, such as premature beats or atrial fibrillation, the accuracy of the readings may be affected.
- Keep the device away from any great power appliances (e.g. high-voltage equipment, X-ray machine, and ultrasound equipment, etc.) to prevent interference during use.
- Keep the device away from active HF surgical equipment, MRI system and RF shielded rooms which generate strong electromagnetic interference.
- Keep the device away from strong electromagnetic fields (e.g. near the mobile power supply, microwave oven, etc.), as it may affect accuracy.
- Do not submerge the device in water or any liquids.
- Do not use the device near flammable anesthetic mixtures.
- Do not stack this device with other equipment to prevent operational issues.
- Keep batteries away from fire and out of reach of children and pets.

1.2 Caution

- For best results, practice multiple times to become familiar with the measurement technique.
- Make sure the batteries are installed correctly before use.
- Inspect the device for any damage before use, especially the cuff. Do not use if damaged to avoid skin injury.
- Verify the device is in proper working condition before use.
- Do not inflate the cuff when not properly wrapped around the arm.
- If discomfort occurs during use, such as prolonged inflation of the cuff, stop the measurement immediately by sliding the power switch to "OFF" or pulling out the air tube plug to release the air.
- Inaccurate measurement values may occur due to unconscious mis-operation. Please follow the provided cleaning instructions in section 11 of this manual to minimize Cross-infection risks for the intended reuse of the product.
- Do not clean the device with volatile solvents.
- Avoid dropping the device or subjecting it to strong impacts.
- Use only specified accessories provided by the manufacturer to maintain proper functionality.
- Do not modify this device. Only replace the battery when needed. Other modifications may create safety risks.
- Follow local laws when disposing of parts or the device at the end of its life cycle.

3

2. Working Principle

The Electronic Blood Pressure Monitor adopts the principle of oscillometric method, receives the blood pressure signal in the armband (cuff) through the pressure sensor, and measures the systolic blood pressure and diastolic blood pressure of the human body after processing by the chip. During the process of measuring blood pressure, the blood pressure sensor detects the number of pulsations of arterial blood vessels, and then converts them into the number of pulsations within a minute to achieve pulse rate measurement.


3. Intended Use

This device is an electronic blood pressure monitor intended for use in measuring blood pressure and pulse rate in adult patient population with the arm circumference range between 8.6 inches to 16.5 inches (22 cm to 42 cm).

4. Product Introduction


Product Name: Electronic Blood Pressure Monitor
Model: BP-205
This equipment is composed of a main body and a cuff. The main body is composed of a central processing unit, a pressure sensor, an air pump, a solenoid valve, a uniform speed vent valve, a PCB board, and an LCD liquid crystal display.

5. Product Structures

4

LCD Display Explanation

6. Symbol Explanation

Graphic Symbol	Content
	WARNING Indicates a potentially hazardous situation, which if not avoided, could result in death or serious injury.
	CAUTION Indicates a potentially hazardous situation, which if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property.
	Please read the instructions before use.
	It belongs to the electrical shock protection BF application part.
IP22	Cover Protection Class IP22: 2 Protected against solid foreign objects of 12.5 mm Ø and greater; 2 Protected against vertically falling water drops when enclosure tilted up to 15°.

5

	The WEEE mark for electronic product recovery indicates that when the end user intends to discard the product, it must be taken to the appropriate facility for recovery and recycling.
	Manufacture date
	Manufacturer
	Production batch number
	CE marking, Notified body number: 0344

7. Product Specification

Item	Specification
Operating Environment	Temperature: 5°C~40°C Humidity: 15%~85% RH Atmospheric pressure range: 70kPa~106kPa
Storage and Transportation Environment	Temperature: -20°C~50°C Humidity: 15%~85% RH Atmospheric pressure range: 70kPa~106kPa
Measuring Method	Oscillometric method
Pressure Sensor	Resistance type pressure sensor
Maximum allowed pressure	300mmHg
Cuff deflation pressure	<10s deflate 260mmHg to less than 15mmHg
Static leakage	≤ 6mmHg for 6min
Manometer test mode	Dynamic pressure mode, Static pressure mode
Display	LCD display
Measuring Range	Blood Pressure/Indication Range: 0-299mmHg (0-39.9kPa) Pulse Rate: 40-199 beats/minute
Measuring Accuracy	Blood Pressure: ±3mmHg (±0.4kPa) Pulse Rate: ±5 %

6

Reproducibility of Blood Pressure Determination	±3mmHg (Root mean square between 10 measurements and their average)
Stability of Blood Pressure Determination	±3mmHg (Maximum deviation between the maximum value of 10 measurements and its average)
Resolution	1mmHg (0.1kPa)
Power Supply	DC 6V (4*1.5V AA)
Memory Capacity	2 users with 99 memories each
Measure Time	<60s
Auto off	Auto power off after 60 seconds of inactivity
Cuff Specification	Approx. 22 to 42cm (8.6" to 16.5")
Dimension	Approx. 106mm (l) × 137mm (w) × 80.9mm (h) (4.2" × 5.4" × 3.2")
Weight	Approx. 290g (batteries and cuff not included)
Package List	Cuff, 4 × batteries, USB-C cable, instruction manual
Lifetime	5 years

Statement

The electronic blood pressure monitor has been performed clinical investigation according to ISO 81060-2: 2018 Non-invasive sphygmomanometers - Part 2: Clinical investigation of intermittent automated measurement type. The mean error and the maximum permissible standard deviation met the requirements of ISO 81060-2: 2018.
-This monitor is clinically investigated in accordance with the requirements of IEC 80601-2-30:2018.
-This monitor's biocompatibility has been verified by test according to ISO 10993 series standards.
-This monitor has been designed and tested in accordance with EN 60601-1:2006+A1:2013 to ensure safety and essential performance.
-Any serious incidents should be reported to the EU representative, manufacturer, and competent authority.

7

8. Settings

8.1 Unit Switching

While the device is on, press the "SET" button to switch units. If the unit is mmHg, press the "SET" button to switch to kPa. When the 0.0kPa is displayed, the setting is successful. Press the "START/STOP" button to save the setting and turn off the device.

8.2 User Switching

While the device is turned off, press the "SET" button twice to enter user selection mode. Press the "SET" button to switch between "user 1" and "user 2".

8.3 Date/Time Setting

While the device is turned off, press the "SET" button 3 times to enter date adjustment mode. The year will start flashing, press the "SET" button to adjust the year.
Press the "SET" button again, the month will start flashing, press the "SET" button to adjust the month.
Press the "SET" button again, the day will start flashing, press the "SET" button to adjust the day.
Press the "SET" button again, the hour will start flashing, press the "SET" button to adjust the hour.
Press the "SET" button again, the minute will start flashing, press the "SET" button to adjust the minute.
Press the "START/STOP" button to save the setting and turn off the device.

8.4 Read the Memory

While the device is turned off, press the "SET" button to view the stored memory values.

8.5 Clear Memory Data

While the device is turned off, press the "SET" button 9 times, the sign "CL" will be displayed. Press the "SET" button and "START/STOP" button at the same time along the cuff with your feet on the floor and arms on the table. When the 0.0kPa is displayed, the setting is successful. Press the "START/STOP" button to save the setting and turn off the device.

Note: All measurement data for all users will be deleted at once.

9. Blood Pressure Measurement Method

9.1 Check and Prepare Before Measurement

- (1) Avoid exercising, bathing, drinking stimulatory drinks (such as coffee or alcohol), and smoking within 30 minutes before the measurement.
- (2) Rest in a comfortable and stable environment for at least 5 minutes before the measurement.
- (3) Remove tight-fitting or thick clothing from your arm, bare the upper arm or wear thin shirts for the measurement.
- (4) Sit in a chair with your feet flat on the floor and arms on the table.
- (5) Make sure the batteries are correctly installed.

9.2 Apply the Cuff on the Arm

- Take out the cuff and open it into a tubular shape, ensuring it is suitable for the arm to fit in.
- Wrap the cuff around the upper arm. Make sure the air tube is positioned along the outer side of the arm. When the cuff is correctly positioned, the bottom edge of the cuff should be 2cm to 3cm (about 1 inch) above the elbow joint.

- Wrap the cuff securely with the Velcro fastener. Leave enough space to fit one finger between the cuff and your arm.

- Plug the arm cuff into the monitor by fully inserting the air tube plug into the air tube port (located on the left side of the monitor).

- Correct measurement posture:
- Rest your elbows on the table.
 - Align the cuff with your heart height.
 - Relax your body with palms facing up.

9

9.3 Start Blood Pressure Measurement

- Press the "START" button while the device is in shutdown or standby mode to initiate the measurement process. The cuff will automatically inflate.

Note: During inflation, please keep your arm still to ensure accurate readings. Any movement may affect the measurement. Feeling the cuff pressure is normal during inflation. Refer to the "Error and Fault Troubleshooting" section for error indications.

- The device will stop automatically after full inflation and proceed to detect your pulse rate.

Once the measurement is complete, your blood pressure and pulse values will be displayed.

Note: Maintain a correct posture and a quiet environment during measurement. To end the measurement, press the "START" button to stop the pressure and deflate the cuff.

- After the measurement, remove the cuff, press the "START" button to manually turn off the device. If no action is taken, the device will turn off automatically after 60 seconds.

Notice:

- Various factors such as measurement site, patient position (standing, sitting, lying down), cuff placement, exercise, and physiological condition can affect blood pressure readings.
- Extreme temperature, humidity, and altitude can impact the blood pressure monitor's performance. Use the device in the environment in accordance with the instruction manual.
- Do not use force or move the arm suddenly during measurement, re-measure.
- Avoid compressing or restricting the air tube.
- Proper cuff tightness is crucial for accurate readings.

10

- Maintain consistency in measurement timing, body position, posture, and device usage.
- Please measure both arms for first-time use.
- If cuff pressure causes discomfort, stop the measurement immediately by sliding the power switch to "OFF" or pulling out the air tube plug to release the air.
- Wait at least 2-3 minutes between measurements to allow the artery to return to the normal state.
- Avoid frequent measurements to prevent blood flow interference and potential injury.
- Kinking of the air tube can lead to continuous cuff inflation, causing harm to the patient.
- Cuff pressurization may temporarily affect the function of other monitoring devices on the same arm.
- Make sure that the blood pressure monitor operation does not excessively impede blood circulation.

10. Error and Fault Troubleshooting

Symbol	Cause	Solution
ER1	Abnormal sensor signal	Check whether the cuff is normal and whether it is worn correctly, re-measure.
ER2	No measurement results	Check whether the cuff is normal, whether it is worn correctly, whether the blood pressure monitor communicates with heart. Using correct sitting posture, keep relax and quiet, re-measure.
ER3	Abnormal measurement results	Check whether the cuff is normal, whether it is worn correctly, whether the blood pressure monitor communicates with heart. Using correct sitting posture, keep relax and quiet, re-measure.
ER4	The cuff is too loose or leaks air.	Check whether the cuff is tight enough to fit one or two fingers into your arm; check whether the air tube joint is inserted into the device; check whether there is obvious air leakage between the air tube and the cuff, re-measure.
ER5	The cuff is too tight or the airway is blocked.	Check whether the cuff is tight enough to fit one or two fingers into your arm; check whether the air tube connected to the cuff is bent or pressed by heavy objects, re-measure.

11

Immunity Test	Compliance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m

Radiated RF IEC 61000-4-3	Compliance
	80MHz- 2.7GHz 10 V/m
	380MHz- 390MHz 27 V/m
	430MHz- 470MHz 28 V/m
	704MHz- 787MHz 9 V/m
	800MHz- 960MHz 28 V/m
	1.7GHz- 1.99GHz 28 V/m
	2.4GHz- 2.57GHz 28 V/m
	5.1GHz- 5.8GHz 9 V/m

Table 3: Not Applicable

Harmonic Emissions (IEC 61000-3-2), Voltage Flicker Emissions (IEC 61000-3-3), Electrical Fast Transients (IEC 61000-4-4), Surge (IEC 61000-4-5), Voltage dips (IEC 61000-4-11), Conducted Immunity (IEC 61000-4-6)	
---	--

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

WARNINGS!

- This product or its parts should not be used in the vicinity of or on the top of other electronic equipment such as cell phone, transmitter or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories sold by power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

14

15. Manufacturer Information


Manufacturer: Shenzhen Brav Electronic Technologies Co., Ltd.
Manufacturer Address: 4-5/F, Building 11, Tongfuyu Industrial Park, Lezhujiao, Huangnabai Community, Hangcheng Street, Bao'an District, Shenzhen, China

CE REP: MedNet EC-REP GmbH
Geschäftsadresse: 10485 Berlin, Germany
E-Mail: service@runstar.com
Telephone: 0162 7868 688

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UK REP: Runstar Healthlink Ltd
85 Great Portland Street/first Floor, London,
United Kingdom W1W 7JT
E-Mail: service@runstar.com

Visit official website or scan QR code to activate lifetime warranty :
www.runstar.com

14. EMC Information and Statement

- ⚠ Precautions:
- This product conforms to the electromagnetic compatibility requirements of IEC 60601-1-2:2014 standard.
- The user shall use the device according to EMC information provided in this manual.
- Portable and mobile RF communication equipment may affect the performance of the product. Avoid strong electromagnetic interference, such as close to mobile phones and microwave ovens.
- The guidance and manufacturer's statement are detailed in the table below.

⚠ Warnings:

This product should not be used close to or stacked with other devices. If it must be used close to or stacked with other devices, it should be observed to verify that it can operate normally under the configuration used.

Table 1: Electromagnetic Emissions

When the 0.0kPa displayed, the setting is success. Press the "START/STOP" button to save the setting and turn off the device.

8.2. User Switching

While the device is turned off, press the "SET" button twice to enter user selection mode.

Press the "" button to switch between "user 1" and "user 2".

8.3. Date/Time Setting

While the device is turned off, press the "SET" button 3 times to enter date adjustment mode. The year will start flashing, press the "" button to adjust the year.

Press the "SET" button again, the month will start flashing, press the "" button to adjust the month.

Press the "SET" button again, the day will start flashing, press the "" button to adjust the day.

Press the "SET" button again, the hour will start flashing, press the "" button to adjust the hour.

Press the "SET" button again, the minute will start flashing, press the "" button to adjust the minute.

Press the "START/STOP" button to save the setting and turn off the device.

13

Table 2: Electromagnetic Immunity

Immunity Test	Compliance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m
	80MHz- 2.7GHz 10 V/m
	380MHz- 390MHz 27 V/m
	430MHz- 470MHz 28 V/m
	704MHz- 787MHz 9 V/m
	800MHz- 960MHz 28 V/m
	1.7GHz- 1.99GHz 28 V/m
	2.4GHz- 2.57GHz 28 V/m
	5.1GHz- 5.8GHz 9 V/m

Table 3: Not Applicable

Harmonic Emissions (IEC 61000-3-2), Voltage Flicker Emissions (IEC 61000-3-3), Electrical Fast Transients (IEC 61000-4-4), Surge (IEC 61000-4-5), Voltage dips (IEC 61000-4-11), Conducted Immunity (IEC 61000-4-6)	
---	--

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

WARNINGS!

- This product or its parts should not be used in the vicinity of or on the top of other electronic equipment such as cell phone, transmitter or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

14

15. Manufacturer Information


Manufacturer: Shenzhen Brav Electronic Technologies Co., Ltd.
Manufacturer Address: 4-5/F, Building 11, Tongfuyu Industrial Park, Lezhujiao, Huangnabai Community, Hangcheng Street, Bao'an District, Shenzhen, China

CE REP: MedNet EC-REP GmbH
Geschäftsadresse: 10485 Berlin, Germany
E-Mail: service@runstar.com
Telephone: 0162 7868 688

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UK REP: Runstar Healthlink Ltd
85 Great Portland Street/first Floor, London,
United Kingdom W1W 7JT
E-Mail: service@runstar.com

Visit official website or scan QR code to activate lifetime warranty :
www.runstar.com

14. EMC Information and Statement

- ⚠ Precautions:
- This product conforms to the electromagnetic compatibility requirements of IEC 60601-1-2:2014 standard.
- The user shall use the device according to EMC information provided in this manual.
- Portable and mobile RF communication equipment may affect the performance of the product. Avoid strong electromagnetic interference, such as close to mobile phones and microwave ovens.
- The guidance and manufacturer's statement are detailed in the table below.

⚠ Warnings:

This product should not be used close to or stacked with other devices. If it must be used close to or stacked with other devices, it should be observed to verify that it can operate normally under the configuration used.

Table 1: Electromagnetic Emissions

12

13. After-sale Service

Visit official website or scan QR code to activate lifetime warranty :
www.runstar.com

A square QR code with a black and white pixelated pattern. It features three distinct square markers in the corners (top-left, top-right, and bottom-left) used for alignment and orientation.

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Aviso:
-Various factors such as measurement site, patient position (standing, sitting, lying down), cuff placement, exercise, and physiological condition can affect blood pressure readings.
-Extreme temperature, humidity, and altitude can impact the blood pressure monitor's performance. Use the device in the environment in accordance with the Instruction Manual.
-Avoid compressing or restricting the air tube.
-Proper cuff tightness is crucial for accurate readings.
-Maintain consistency in measurement timing, body position, posture, and device usage.

- Please measure both arms for first-time use.
- If cuff pressure causes discomfort, stop the measurement immediately by sliding the power switch to "OFF" or pulling out the air tube plug to release the air.
- Wait at least 2-3 minutes between measurements to allow the artery to return to its normal state. Frequent measurements can harm the patient.
- Avoid frequent measurements to prevent blood flow interference and potential injury.
- Kinking of the air tube can lead to continuous cuff inflation, causing harm to the patient.
- Cuff pressurization may temporarily affect the function of other monitoring devices on the same arm.
- Make sure that the blood pressure monitor operation does not excessively impede blood circulation.

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Fabricante: Shenzhen Brav Electronic Technologies Co., Ltd.
EC REP MedNet EC-REP GmbH
EU REP Runstar GmbH
UK REP Runstar Healthlink Ltd