VITRO

ECG Smart Band

User Manual

Navigation Page

English - - - - - (1~14)

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Español - - - - - (41~54)

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
Product Introduction

VITRO is a premium screenless smartwatch equipped with ECG functionality, integrating health monitoring, fitness tracking, and smart reminders to provide users with comprehensive health management.

How to Use

Vibration reminder

- **Power On / Power Off:** Long press the side button to power on, Vibrates once (250 ms)

Pairing with Phone


1. Download and install the "G Band" App (available for iOS and Android).
2. Turn on Bluetooth on your phone and search for the watch via the App.
3. Follow the instructions in the App to complete pairing.

Charging


1. Use the original magnetic charging cable to charge the watch.
2. Avoid using non-original chargers to prevent device damage.

Wearing

Wear the watch 3–2 cm above your wrist bone, ensuring the sensors maintain firm contact with your skin for accurate measurements.

ECG Detection

Before testing blood pressure and heart rate, you need to rest quietly for more than 5 minutes. Do not eat, smoke, take medicine or consume caffeine, etc. 2 hours before the test.

During the test, adopt a comfortable and quiet sitting posture, maintain a calm state of mind, and refrain from speaking.

When testing, please wear the wristband on your left wrist. Then, gently press the button at the bottom of the wristband with your right middle finger. During the test, the two electrode contacts at the bottom should be in contact simultaneously. Do not press the electrode sheet too hard, as shown in the following picture: If the test fails, please breathe as per the instructions or moisten the skin at the contact point with salt water and then retest!

Vibration Reminder Indicators

1. Call Reminder Setting

- The call reminder function must be enabled in the App; otherwise, the watch will not vibrate for incoming calls.
- By default, the call reminder feature is enabled in the App.

2.Call Reminder Rule


- When an incoming call is detected, the watch vibrates for 250 ms with a 750 ms pause in a continuous loop until the call is answered, ended, or the watch button is pressed.

3. Health & Message Reminder Rule

- The watch vibrates four times consecutively to notify health alerts or incoming messages.

Features Overview

1. Daily Activity Steps Summary

- Tracks total daily steps (walking, running, fitness walking).
- Default target: 8000 steps, adjustable via the app.

Calories Summary

Records total calories burned throughout the day.

Distance Summary

Tracks total movement distance in 24 hours.

2. Health Monitoring


Sleep Monitoring

The device automatically detects and analyzes your sleep quality and duration, including long sleep and short naps.

- Long Sleep: Records deep sleep, light sleep, REM, awake periods, sleep onset time, and wake-up time.
- Short Naps: Records nap duration and time periods.
- Total sleep duration and daily comparison are available in the App. You can also enable "Scientific Sleep Mode" in the App for more detailed sleep analysis reports.

Heart Rate Monitoring

24-hour heart rate tracking with abnormal heart rate vibration alerts.

- Can be toggled on/off in App.
- Measurement interval: 1 min (during sleep) / 10 min (normal); uploads every 30 min.

High Heart Rate Alert

When the heart rate exceeds the set threshold, the watch vibrates 5 times as a warning.

Blood Pressure

- Optional automatic monitoring (every 10 min).
- Supports manual calibration mode for higher accuracy.


Blood Oxygen (SpO₂)

Monitors oxygen saturation levels for reference.

- Same measurement frequency as heart rate.
- Data uploaded every 30 minutes.

HRV | Heart Rate Variability

Analyzes physical recovery and stress resilience.

Range:

- 1–29: Low
- 30–60: Normal
- 61–101: Good
- 101+: Excellent


Blood Glucose

24-hour monitoring with measurements every 10 minutes.

- Calibration is manual (optional).

Stress Monitoring

Monitors stress levels in real-time.

Range:

- 1–29: Relaxed
- 30–59: Normal
- 60–79: Moderate
- 80–99: High

ECG Monitoring

Performs real-time ECG analysis and heart rhythm detection.

Risk Level Range:

- 1–5: No Risk
- 5–50: Low Risk
- 50–70: Medium Risk
- 70–99: High Risk

Body Composition

BMI Classification:


- <18.5: Underweight
- 18.5–24.0: Normal
- 24.0–28.0: Overweight
- ≥28.0: Obese

Menstrual Cycle Tracking

Set your cycle in the app; the watch displays the corresponding status and reminders.


Exercise Features

VITRO supports 188 exercise modes, covering a wide range of activities:

- Running, cycling, swimming, yoga, dance, ball games, and winter sports.
- Supports both indoor and outdoor workouts.
- Exercise data (steps, heart rate, calories, routes) can be viewed in real time via the app.


SOS & Fall Detection

(After setting up the SOS emergency contact, you need to light up the SOS icon to turn it red.)

1. Certain models include SOS emergency and fall detection functions.
2. When a fall is detected or SOS is triggered, the watch will:
3. Automatically call your emergency contact.
4. Vibrate continuously to notify you.

Relatives and Friends Mode

To use the Relatives and Friends Mode you need to enable Data Cloud Sync on the Data Services page first, then enter the mode and add member accounts via QR code scanning or email address.

Smart Reminders

-  Incoming call and message notifications (via vibration)
-  Sedentary alerts
-  Alarm reminders
-  Daily goal achievements
-  Health alerts (e.g., abnormal heart rate or low oxygen)

FAQ

1. Watch Won't Power On

- Ensure sufficient battery and press the side button for 5 seconds.

2. Watch Failed to Pair with Phone

- Check that Bluetooth is enabled and the app is up-to-date.

3. ECG Reading Abnormal

- Keep the watch firmly on your wrist and keep it still; repeat the measurement if needed.

Safety & Maintenance

Safety

1. Do not use in high temperature, high humidity, or strong magnetic environments.
2. If you experience any skin abnormalities or allergies, stop wearing the watch and

consult a doctor.

3. Do not use in high temperature, high humidity, or strong magnetic environments.
4. Stop wearing if skin irritation occurs and consult a doctor.

Maintenance

1. Use a soft cloth to clean the watch face and sensor; avoid using chemical solvents.
2. Avoid dropping, collision, or prolonged water exposure (not for diving).
3. Avoid prolonged exposure to sunlight to prevent aging of the strap and case.
4. Clean with a soft cloth; avoid chemical solvents.
5. Avoid dropping, collision, or prolonged water exposure (not for diving).
6. Avoid prolonged sun exposure to prevent aging of the strap and case.

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FCC Warning:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 0cm between the radiator and your body.