

产品名称：**M1磁吸款光疗面罩说明书**

产品尺寸：**160 * 140 mm**

产品材质：**157g 双铜版，骑马订**
封面封底过哑膜

产品工艺：**四色**

设计师：**曾凯欣**

日期：**20250721**

备注：

运营：**陈媛媛**

产品：**冯海烽**

版本：**V1**

本品牌及其所使用的标志、名称、图像以及所有内容设计，及其
版权均均为本公司所有，未经授权，任何组织、个人、商业团体不
得擅自使用、翻印，本公司保留追究其法律责任的权利！





LED LIGHT THERAPY MASK



FOREWORD	1
CONTROLLER	2
WEARING METHOD	3
INTRODUCTION TO ILLUMINATION	5
Instructions for use and treatment recommendations:	5
PRODUCT PARAMETERS	8

FOREWORD

Thank you for choosing our products!

Simplify your life with LED Face Photon skin rejuvenation!


Please read the following instructions carefully before using the product, as they may explain important operating procedures.

LED Facial Beauty Device photon skin rejuvenation is a portable and rechargeable LED Facial Beauty Device photon beauty device for home use. The surface structure of this product is made of food-grade silicone material, which is convenient and easy to use, foldable and portable, and can be fixed with adjustable nylon buckle, suitable for different groups of people.

It always has face 88x4=352 rays, It basically solves the problem of whitening and skin ageing, that is: increase cell activity; promote cell metabolism, so that the skin secretes collagen and fibrous issue, to make it self-accelerate blood flow back, increase skin elasticity, improve skin after, all conditions, so as to achieve whitening, skin adjustment, anti-wrinkle, light spot, and other repair effects.




CONTROLLER



Turn on and off
Long press for 2 seconds



Switch Color
Touch to click



Switch usage time
Click twice



Adjust brightness
Touch to click

PS: Please connect the control unit to the mask before operation.


Photosensitivity test

Before using this product for the first time, please carry out a light sensitivity test:

1. Place the device on your arm, switch it on and let it glow for 10 minutes. If symptoms such as redness or itching occur, you may have photosensitive skin and it is recommended that you stop using the device or consult a beautician before using it again.

WEARING METHOD

When wearing the mask, please note the position of LEFT and RIGHT as shown in the picture.



2. Please cleanse your face before use and keep it dry (without make-up) to avoid photosensitive ingredients that may be contained in cosmetics and impair the effectiveness of skin care.

3. Reddening of the skin may occur after application, which is a normal skin reaction and will normalise within 1-2 hours.

If you experience a burning sensation or other discomfort during application, discontinue use immediately.

Charging instructions:

1. This product contains lithium-ion polymer batteries and should not be used while charging (Note: If used while charging, the battery and device will be damaged and the device will not light up and function)

2. Before first use. Please fully charge the device. When charging, connect the USB charging cable to the control unit, insert the plug, switch on the device and start charging.

3. When the indicator light flashes, the device is charging; when the indicator light is off, charging is complete.

Tips:

Before use, please irradiate the face for 1 minute, no discomfort symptoms can continue to use.

Please start with a low gear, build up the skin tolerance and choose the energy level according to the actual skin condition. Please do not blindly increase the time and frequency of light, you need to adjust the length of light according to the actual condition of the skin. It is recommended to use this product 3 - 5 times a week, the one-time application time is 5/10 minutes. Once results are visible, reduce the application time to 2-3 times a week.

WARNING

- Do not use the device while it is charging.
- Please do not fold up the mask measuring device for storage, otherwise the cables inside will be damaged.

INTRODUCTION TO ILLUMINATION

Red light (630nm):

It improves cell activity, accelerates metabolism and promotes collagen formation, lightens spots, rejuvenates skin and eliminates wrinkles, increases skin elasticity, improves skin condition, repairs damaged skin, reduces pores and tightens skin, has anti-ageing and antioxidant effects.

Instructions for use and treatment recommendations:



1. Photosensitivity test (must be carried out on first application)

Select an inconspicuous area (e.g. the wrist or inside of the arm), switch on the dim light of the mask (level 1) and irradiate for 5-10 minutes.


Within 24 hours, watch for abnormal reactions such as redness, swelling, itching, etc. If you feel unwell, please discontinue use and consult a doctor.



2. Cleanse the skin

Wash your face and hands thoroughly with warm water before use to ensure that no make-up, oil or dirt remains.

Dry the skin carefully with a towel to keep it dry (residual moisture can impair the effect of phototherapy).



3. Connect the controller to the mask



4. Select the light color:

Choose the appropriate light mode according to your skin care needs.



5. Wear the mask

Adjust the elastic band so that the mask is completely against the face to prevent light from escaping.

Connect the control unit and check that the indicator light lights up normally.

6. Recommendations for frequency and duration of use

This mask has 3 intensity levels. Please start with the lowest level and adjust it gradually.



First week (adaptation phase)

Level: 1 (it is recommended to start with level 1)

Frequency: 3 times per week

Duration of a single session: 10 minutes

Note: It is normal to feel a slight fever during the first application. If you feel a tingling sensation or severe discomfort, stop immediately.



Second week (strengthening phase)

Level: 2 (adjusted according to skin tolerance)

Frequency: 3 times per week

Duration of a single session: 10 minutes (extended to 20 minutes if no discomfort occurs)

Note: For sensitive skin, a low level can be maintained and long-term exposure avoided.



Third and fourth week (consolidation phase)

Level: 3 (flexibly adjusted depending on skin condition)

Frequency: 3 times per week

Duration of a single session: 10 minutes

Note: For dry or reddened skin, the duration can be shortened or applied at intervals.



7. Care after application

Immediate reaction:

Skin may be temporarily reddened or dullness (disappears within 1-3 hours) and may take 1-3 day to recover after blue light exposure.

There may be a feeling of tightness after the first application, which is a normal phototherapy effect.



8. Skin care instructions:


Apply moisturising or repairing skin care products (e.g. essence, cream) after application to enhance the effect.

Avoid the immediate use of irritating ingredients (e.g. fruit acid, retinol).

PRODUCT PARAMETERS

Name	LED Light Therapy Mask
Number of lamp beads	352 high-quality LEDs
Charging Voltage	5V 1A

Charging time	About 100 min
Net weight	150g
Battery capacity	1500mAh



- Beauty instrument x 1
- Fixing strap x 2

- Controller x 1
- Charging cable x 1

- Eye protection x 2
- Manual x1

FCC Caution.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

* RF warning for Portable device:

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.



V1 **FC** **CE**  
Made in China