

300005273225 (13/8/2026)

Philips Oral Healthcare, LLC.
22100 Bothell Everett Highway
Bothell, WA 98021
U.S.A.
1-888-744-5477
www.philips.com/Sonicare
www.philips.ca/sonicare

©2026 Koninklijke Philips N.V.
(KPNV). All rights reserved.
Philips and the Philips shield are
trademarks of KPNV.
Sonicare and the Sonicare logo
are trademarks
of Philips Oral Healthcare, LLC
and/or KPNV.

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English

Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Important safety information

Only use the product for its intended purpose. Read this information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries. Accessories supplied may vary for different products.

WARNINGS

- This appliance can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
- Children should be supervised to ensure that they do not play with the product.
- Do not clean any part of the product in the dishwasher unless otherwise instructed in the Cleaning section.
- This appliance contains no user-serviceable parts. If the appliance is damaged, stop using it and contact the Consumer Care Center in your country.
- This appliance has only been designed for cleaning teeth and gums.
- Replace the brush head every 3 months or sooner if signs of wear appear. Stop using a brush head with crushed or bent bristles.
- Rinse the handle, especially the brush head connection. Gently clean around the rubber seal.
- Make sure brush head and toothbrush are dry before storing in travel case.
- Avoid direct contact with products that contain essential oils or coconut oil. Contact may result in bristles dislodging.

MEDICAL WARNINGS

- Stop using the appliance and consult your dentist/doctor if excessive or prolonged bleeding occurs after use, or if you experience discomfort or pain during use.
- If you have a pacemaker or other implanted device, contact your doctor or the manufacturer of the implanted device prior to use.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.

Battery safety instructions

- 1 X AAA battery required.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, or changes color, stop using the product and contact Philips.
- Do not place products and their batteries in microwave ovens or on induction cookers.

- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.
- If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.
- When you handle batteries, make sure that your hands, the product and the batteries are dry.
- To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal objects (e.g. coins, hairpins, rings). Do not wrap batteries in aluminum foil. Tape battery terminals or put batteries in a plastic bag before you discard them.
- Only operate this product with type of disposable battery specified in the user manual.
- Always keep AAA batteries out of reach of children. Swallowing these batteries can result in serious injury (chemical burns, perforation of soft tissue) or death. Severe burns can occur within 2 hours of ingestion. If swallowed, seek medical attention immediately.
- Do not recharge non-rechargeable batteries.
- Insert batteries with the + and - poles pointing in the direction indicated in the battery compartment or holder.
- Remove batteries from the product if you are not going to use it for some time.
- Do not leave empty disposable batteries inside the product.
- Do not mix different brands and types of batteries, do not mix new and used batteries and do not use batteries with different date codes.
- Do not use non-rechargeable batteries in place of rechargeable batteries.
- Do not use modified or damaged batteries.

Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

SAVE THESE INSTRUCTIONS

Explanation of symbols

The following symbol may appear on the product:
Read Operator's Manual (Fig. 2): This symbol means read the instructions for use, before using the appliance.

Intended use

Philips One power toothbrushes are intended to remove adherent plaque and food debris from the teeth to reduce tooth decay and improve and maintain oral health. Philips One power toothbrushes are intended for consumer home use. Use by children should be with adult supervision.

Your Philips One (Fig. 1)

- 1 Travel case door
- 2 Travel case
- 3 Brush head
- 4 Handle
- 5 Power button
- 6 Removable end cap with coin slot
- 7 Replaceable AAA-battery

Note: Accessories included may vary based on the model purchased.

Getting started

Attaching the brush head

- 1 Align the brush head so the bristles point in the same direction as the front of the handle. Press the brush head down onto the metal shaft (Fig. 3).

Note: It is normal to see a slight gap between the brush head and the handle (Fig. 4). This allows the brush head to vibrate properly.

Battery replacement

Your power toothbrush is powered by a replaceable AAA-battery. A single AAA battery should provide enough power for two 2-minute brushing sessions every day for 3 months. If used more often, the battery will drain before 3 months.

Note: Make sure your hands and the product are dry when you insert the battery.

- 1 Insert a coin into the coin slot at the bottom of the handle and turn counterclockwise (Fig. 5). Gently shake the battery out of the handle.
- 2 Insert a fresh AAA battery into the battery compartment (Fig. 6). Make sure the + and - poles of the battery point in the right direction.
- 3 Place the removable end cap back on the handle and use a coin to turn clockwise until it is secure and the lock indicator markings on the handle and end cap line up (Fig. 7).

Note: If the indication on the handle and end cap do not line up, it may result in your toothbrush not working

Travel case

Your Philips One comes with a travel case to store your power toothbrush when traveling.

- 1 Place your Philips One toothbrush with the bristles facing downwards in the travel case and snap closed.
- 2 To remove, open and press down on the bottom of the handle (Fig. 8).

Using your Philips One toothbrush

If you are using this toothbrush for the first time, it is normal to feel more vibration than when using a manual toothbrush. The Philips One uses microvibrations to help gently clean your teeth. Please follow the brushing instructions below for the best experience.

Philips One comes with a SmarTimer and QuadPacer to help you brush evenly throughout your mouth.

SmarTimer

The SmarTimer indicates that your 2-minute brushing cycle is complete by automatically turning off the toothbrush at the end of 2 minutes.

Note: Dental professionals recommend brushing no less than 2 minutes twice a day.

QuadPacer

The QuadPacer is an interval timer that beeps every 30 seconds to remind you to brush all teeth evenly by dividing your mouth into 4 sections (Fig. 9). Begin brushing in section 1 (outside upper teeth) and brush for 30 seconds before you move to section 2 (inside upper teeth). Continue brushing in section 3 (outside lower teeth) and brush for 30 seconds before you move to section 4 (inside lower teeth).

Note: If you press the power button after you started the brushing cycle, the toothbrush will pause. After a pause of 30 seconds, the SmarTimer and QuadPacer resets.

Brushing instructions

Please follow the brushing instructions for the best experience.

- 1 Wet the bristles and apply a small amount of toothpaste.
- 2 Place the toothbrush bristles against the teeth at a slight angle (45 degrees) (Fig. 10), pressing gently to make the bristles reach the gumline or slightly beneath the gumline.
- 3 Press the power button to turn brush on.
- 4 Brush your teeth with back and forth strokes so that the bristles reach all surfaces of the teeth. Change to another section when the QuadPacer beeps. Continue this motion throughout the 2-minute brushing cycle.

Note: The bristles should slightly flare.

- 5 To clean the inside surfaces of the front teeth, tilt the brush vertically and make several up-and-down strokes (Fig. 11).
- 6 After you have completed the brushing cycle, you can spend additional time brushing the chewing surfaces of your teeth and areas where staining occurs (Fig. 12). You can also brush your tongue, with the toothbrush turned on or off, as you prefer.

Your power toothbrush should be safe to use on braces (brush heads wear out sooner when used on braces) and dental restorations (fillings, crowns, veneers) if they are properly adhered and not compromised. If a problem occurs, please follow up with a dental professional.

Battery status

This product does not come with a battery status indicator. As the battery gets weaker, you may experience a drop in vibration power. This is normal. While the microvibrations are still effective, you can increase the power by replacing your battery (see 'Battery replacement').

Cleaning

You should clean your toothbrush and accessories regularly to remove toothpaste and other residue. Failure to clean your product may result in an unhygienic product and damage to it may occur.

Warning: Do not clean product or accessories in dishwasher, microwave, with chemicals or with boiling hot water.

Toothbrush handle

- 1 Remove the brush head and rinse the metal shaft area with warm water (Fig. 13). Make sure you remove any residual toothpaste.
- 2 Wipe the entire surface of the handle with a damp cloth.

Note: Do not tap the handle on the sink to remove excess water.

Brush head and travel case

- 1 Rinse the brush head and bristles after each use and allow to air dry.
- 2 Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week. Rinse the travel case weekly with warm water and wipe clean with a damp cloth.

Storage

If you are not going to use the product for a long time, clean it (see 'Cleaning') and store it in a cool and dry place away from direct sunlight.

Warning: If product will not be used on a regular basis, remove the battery to prevent possible damage due to chemical leakage. If battery does leak, remove it carefully. Do not allow bare skin to touch leaking fluid.

Replacement

Brush head

Replace the brush heads every 3 months to achieve optimal results. Only Philips One for Kids brush heads are compatible with the Philips One for Kids toothbrush.

Locating the model number

Remove brush head from handle. Look at the top of the handle near the metal shaft for the model number (HY110X).

Troubleshooting

This chapter summarizes the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Possible Cause	Solution
I cannot fully attach the brush head. There is a gap between the brush head and handle.		This gap is normal and is needed for the brush head to vibrate appropriately. The brush head must be able to move to create the right amount of vibrations.
I cannot fully attach the brush head at all.	A standard Sonicare brush head is being used.	Standard Sonicare click-on brush heads will not fit on the toothbrush. Only use brush heads specifically created for the Philips One with a model number that begins with BH1022.
My Philips Sonicare toothbrush vibration is less powerful than before.	The brush head is too close to the body of the toothbrush or the toothbrush needs a new battery.	Remove the brush head from the handle and re-set it on the metal shaft leaving a small gap. (1-2 mm). Ensure a new AAA battery is inserted following the battery replacement instructions under 'Getting started'.

Full Two-Year Warranty

Philips warrants its products for a specific period of time after the first date of purchase. If the product is purchased in the USA, Philips North America LLC (USA) or, if the product is purchased in Canada, Philips Electronics Ltd., warrants this new product against defects in materials or workmanship for a period of two years from the date of purchase (as evidenced by a valid proof of purchase), and agrees to replace any defective product without charge. Philips will refund the purchase price if it is unable to provide replacement.

Parts that are subject to normal wear and tear, accessories, replacement batteries and consumables such as toothbrush heads, are not covered by the terms of this warranty.

IMPORTANT: This warranty does not cover damage resulting from accident, misuse or abuse, lack of reasonable care, the affixing of any attachment not provided with the product or loss of parts or subjecting the product to any but the specified voltage.*

NO RESPONSIBILITY IS ASSUMED FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES.

For more information and in order to obtain warranty service, simply go to www.philips.com/support. NOTE: This warranty gives you specific legal rights, and you may also have other rights which vary from state to state, province to province or country to country. Some states do not

allow the exclusion or limitation of incidental or consequential damages, so the above exclusion and limitations may not apply to you. This warranty is transferable. *Read enclosed instructions carefully.

THE WARRANTY DOES NOT COVER THE FOLLOWING:

- Indirect or consequential damages, such as loss of time, loss of data, or loss of income. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.
 - Compensation for activities done by yourself, such as regular maintenance, installing firmware updates, or saving data.
- ## THE WARRANTY DOES NOT APPLY:
- If the proof of purchase, model number, serial number, or production date code on the product has been made illegible.
 - If the product was purchased as faulty stock.
 - If the defect is the result of abuse or excessive use outside the intended purpose.
 - If the defect is caused by environmental conditions, connecting peripherals, additional equipment, or accessories (including software) that are not in accordance with the user manual.
 - If the unit has been damaged, including but not limited to damage by animals, lightning, abnormal voltage, fire, natural disaster, transportation, or water (unless the user manual expressly states that the product may be rinsed).
 - For non-electrical parts, normal wear and tear, accessories, replacement batteries or replacement of consumable parts such as toothbrush heads.
 - If the product does not function properly because it was not originally designed for use in the country where the product is used, which may occur if you import the product.
 - If the product does not function properly due to problems with connection with service providers, such as for the internet.

Disposal

Your product is designed and manufactured with high quality materials and components, which can be recycled and reused. For recycling information, please contact your local waste management facilities or visit: www.recycle.philips.com.

- This product contains batteries:
- Dispose of batteries properly. Do not incinerate. Batteries may explode if overheated.
 - Do not wrap in metal or aluminum foil. Tape the waste battery terminals before discarding.
 - It is suggested that you contact your local town or city to determine proper battery redemption site(s) in your area.
 - This symbol means that electrical products and batteries shall not be disposed of with normal household waste (Fig. 14).
 - Follow your country's rules for the separate collection of electrical products and batteries.
 - To remove the battery, follow the instructions in the battery replacement chapter.

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