
Cordless Aesthetic

FACE SONO

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Manual

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Homecare Device

Thank you for purchasing a product from JLU Co., Ltd.

This product is a beauty device, not a medical device.

Please note that the user experience and perceived effects may vary by individual. If any abnormal symptoms occur, discontinue use immediately

Table of Contents

01. Caution	-----	13
02. Product Introduction		
- Product Components & Name	-----	14
- Product Specifications & Charging Information	-----	15
03. Usage Method		
- Usage Method	-----	17
04. Warranty & A/S		
- Warranty & A/S	-----	20
- Product Warranty & Product Registration Guide	-----	21


※ CAUTION

1. This product is a beauty device, not a medical device. Please refrain from using it for any purpose other than beauty care.
2. Remove the Protective film from the head before using the device.
3. Before using the product, please make sure it is functioning properly.
4. Avoid use on children·pregnant·breastfeeding women.
5. Please keep out of reach of children.
6. Avoid use on people who require health precautions, such as having pacemaker·device implants·history of severe illness·outpatient treatment.
7. Remove metal accessories such as necklaces·earrings etc., before use.
8. Avoid use on areas with skin diseases such as wounds·infections·boils.
9. It is recommended to consult with a medical professional when using with dental treatments such as implants·nerve treatments·orthodontics.
10. Avoid use on pupils(eyes), around the eyes and where the thyroid gland passes.
11. This product is a skin care device intended for use on the face only. Avoid using it on any area other than the face.
12. Before applying the product to the customer's skin, always check the condition of the contact surface first (temperature, foreign substances, etc.).
13. During use, make sure the device is in full contact with the skin.
Do not leave any space between the device and the skin.
14. Use carefully, as prolonged use on one area may cause skin irritation.
15. Set the appropriate power based on the skin sensitivity.
16. Within 7 days after treatment, avoid activities that may irritate the skin, such as sauna·steaming·exfoliation
17. We are not responsible for any side effects that may occur as a result of not following our recommendations.


※ Product Care

1. If a detachable power supply is provided, only use the device with the power supply provided with it.
2. This product does not include a power supply unit. Use a power supply unit that complies with the product specifications.
3. Do not modify or disassemble·repair the product. If the user modifies·disassembles·repairs the product, free A/S service cannot be provided within the warranty period, therefore be sure to use the designated A/S center.
4. Do not touch with wet hands, do not wash or immerse the product in water. Be careful as it may cause electric shock·malfunction.
5. Do not store the product in humid environments (bathrooms·laundry rooms·saunas etc), as this may cause malfunction. Store it in a dry place, away from moisture.
6. Avoid dropping the product or subjecting it .
7. After use, gently wipe the head(contact area) with alcohol to maintain hygiene. Use a soft cloth for disinfection, as rough cloths or strong wiping can cause scratches or damage.
8. Use wet tissues or a soft cloth to clean the main body, as using alcohol may cause discoloration·contamination.
9. For detailed product care·A/S services, follow the manufacturer's warranty regulations.

Components


Main Body

Lymfit Hydra
Sonic Essence Gel 300g
(Complimentary)

Name

Product Specifications

*This product's specifications may be subject to change without separate notice for performance improvement.


*The specifications of this product may vary depending on the measuring tools, location, and other measurement conditions.

Model Name	JAH-S100A	Product Name	Skin care device
Product Size	H220 x D47 mm	Product Weight	About 370g
Input	DC 5V 2A	Battery	2400mAh
Charging Cable	C-TYPE	LCD	1.14inch (135 x 240 px)
Charging Time	About 180 minutes	Usage Time	About 150 minutes
Manufacturer/Seller	JLU Co., Ltd.	Date of Manufacture	Separate Notation
A/S Center	+82 10-4505-7989	Customer Center	+82 51-331-7185


Charging Information

When connecting the charging cable to the device or cradle, make sure the connector is fully inserted. An incomplete connection may result in improper charging. If the device does not charge or power on, please check the charging connection.


Device Charging

Cradle Charging


Charging Screen

Charging Screen

The icon will blink while charging.

Charging Complete Screen

Charging Information

1. We recommend using an adapter with the rated voltage (DC 5V, 2A).
Charging with a power bank is also possible, but the charging speed may be slightly slower.
2. After a full charge, when the power is turned on, the battery indicator may show 3 bars instead of 4.
3. This does not mean the battery is undercharged. The 4 bar section is designed to indicate only the very highest voltage range, so it may drop to 3 bars after just a short period of use.
4. The battery is still fully charged, and this does not affect the actual usage time.
5. Near full charge or right before complete discharge, the device is designed to charge slowly with a micro-current for battery safety. Because of this, charging may feel slightly longer than expected.
6. If the charger is reconnected after a full charge, the display may still show "CHARGING", which is also normal. You can use the device with confidence.
7. Do not operate the device while charging (e.g., pressing buttons), as this may cause malfunction.
8. If the charging port is not properly aligned, or if dust or foreign substances accumulate on the connector, poor contact may occur. Please check carefully before charging.

Device Care Method

Alcohol Disinfection Allowed

Use a soft cloth for disinfection

Alcohol Prohibited


Using a rough cloth or wiping the surface too hard may cause scratches or damage. Product damage due to user negligence is not eligible for exchange or refund.

*Use wet tissues or a soft cloth to clean the device body, as using alcohol may cause discoloration-contamination.

Usage Cycle	Can be used daily, recommended for 10~15 minutes a day *Depending on your skin condition, adjust the intensity and frequency of use.
Prohibited Areas of Use	Pupils(eyes), around the eyes the area where the thyroid gland passes
Mode	3 types (QUAD-1/3/10/17MHz cross, DUAL-3/10MHz cross, MASSAGE-3MHz, MOIST-10MHz)

Gyro Sensor Function

To prevent malfunction and overheating, the device is equipped with a movement detector gyro sensor. If no movement is detected during operation, the output will automatically stop. Once movement is detected again while the gyro sensor is active, the output will automatically resume.

Precautions

*Precautions for Water·Moisture, and the Use of Ultrasound Gel·Skincare

1. This product is not waterproof. Please avoid storing or using it in humid environments such as bathrooms, laundry rooms, or saunas.
2. Exposure to water or moisture may cause product malfunction.
3. If ultrasound gel or skincare are applied directly onto the device, the contents may seep into gaps in the product or into the charging port, which may cause malfunction. Any damage-loss resulting from this will be considered consumer negligence, and the manufacturer will not be responsible. Apply ultrasound gel or skincare to the face first before using the device.

Usable

Before using the device,
apply ultrasound gel or cosmetics to your face.

* Use cosmetics with a light,
liquid texture rather than cream-type products.

* Use sheet masks instead of gel-type masks.

Prevent liquid from entering the device

Gaps in the product

Vent


Charging port

03. Usage Method - Quick Guide

19

STEP 01

After cleansing, dry your face and apply a sufficient amount of the gel to the area to be treated.

STEP 02

Select the desired mode(frequency) and intensity level, then start.

STEP 03

Rub each area in order 1,2,3 for about 30 seconds each, then gently handle the device in the direction of 3→2→1.

STEP 04

Rub the entire face starting from the chin area then cheeks→ eyes→forehead.


STEP 05

After facial treatment, gently handle the device from the center of the face to 3→2→1.

Operation Method

POWER ON	Press the power button for more than 2 seconds and the device will turn on with the screen.
POWER OFF	When it turns on, press the power button for more than 2 seconds to turn it off.
Mode Settings	Set the frequency by briefly pressing the mode button. (QUAD, DUAL, MASSAGE, MOIST)
Power Settings	Set the power by briefly pressing the power button. (Level 1~5)
Pause	It is possible to pause with the start/stop button.

Screen Information

- ① Battery status
- ② The currently set mode
- ③ The currently set power
- ④ Pause icon
- ⑤ Remaining Time

This warranty policy applies only within the warranty period and conditions specified by our company. For A/S service, please contact the place of purchase or our headquarters A/S center.

※ Warranty Period

Within 1year from the date of manufacture, free maintenance is provided upon product registration. If a manufacturing defect in the main body is confirmed during the warranty period, free A/S service will be provided.

※ Paid A/S

① Products outside the warranty scope

- Products past the warranty period
- Battery performance degradation and wear of consumables due to normal use

② Consumables and accessories

- Consumables (head, cartridge, cable, adapter, handpiece, cart, etc.)
- Separately provided accessories (stand, cradle, etc.)

③ Customer negligence or improper management

- Damage caused by customer negligence such as dropping·impact·water exposure
- Damage due to improper storage or usage conditions
- Exterior damage (scratches, discoloration, peeling of coating, etc.)

④ Unauthorized modification·environmental factors

- Failures caused by unauthorized disassembly·modification·repair
- Issues due to unofficial software·firmware changes
- Damage caused by natural disasters (fire, earthquake, flood, etc.)

※ Exceptions (A/S is not available)

Lost or stolen

※ A/S Shipping Fee

For services conducted via courier, shipping costs are shared between the seller and the customer.

- Before repair : Customer bears the cost
- After repair : Seller bears the cost

※ A/S Center

- Contact Number : +82 10-4505-7989
- Operating : 9:00 AM ~ 6:00 PM Closed on Sat/Sun/public holidays(Lunch break 12:00 PM - 1:00 PM)
- Address : JLU Co., Ltd. 411 office, Busan Knowledge Industry Center, 111, Hyoyeol-ro, Buk-gu, Busan, Republic of Korea

Product Warranty

Thank you for purchasing JLU products.

This product is a genuine product manufactured by the headquarters and shipped after passing strict quality control and inspection processes.

Free service within the warranty period is available if you complete the product registration through the homepage.

A/S Center : +82 10-4505-7989

[How to Register the Product]

Visit the official website (www.jlu.kr) → Go to the Product Registration page
→ Log in → Enter your product serial number

*S/N : Separate notation on the product

FCC Compliance Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates,uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

How to check the electronic label (e-label)

This product utilizes an electronic labeling (e-label) method to display conformity assessment information on the product screen.

How to access the e-label


After turning on the device, please check the e-label information by following the steps below.

- ① Please press and hold the Mode button for 3 seconds.
- ② After 3 seconds, the conformity assessment information is displayed on the E-label screen.

Exit View

Releasing the Mode button returns the display to the operation screen.

※ E-label information can be viewed directly on the product screen without entering a password or connecting additional devices.

※ Certification information for this product is provided not only via the on-screen E-label but also on the packaging box and in the manual, allowing you to check the details at any time.