Smart Body Fat Scale

Model No.: FI2419

User Manual

Read this manual before use and save for future reference

Contents

ENGLISH	1-14
DEUTSCH	15-28
FRANÇAIS	29-42
ITALIANO	43-56
ESPAÑOL	57-70

Visit <https://www.runstar.com> or scan the QR code below to register and get lifetime warranty:

Thank you for purchasing Runstar Smart Body Fat Scale, your ultimate choice for comprehensive body analysis. This advanced scale provides unparalleled convenience by measuring an extensive range of vital health indicators. With the smart body scale, you can effortlessly monitor your body weight, BMI, body fat rate, fat-free body weight, subcutaneous fat, visceral fat, body water, skeletal muscle rate, muscle mass, bone mass, protein levels, BMR, body age, and much more. Experience a new level of precision and insight into your overall well-being with the Runstar.

IMPORTANT: Prior to using the device, please carefully review the instructions in this manual. Failure to follow the instructions or using the scale in ways not specified may lead to personal injury or damage to the scale.

WARNING

- ★ DO NOT use this scale if you have medical implants, such as pacemakers.
- ★ DO NOT step on the edge of the scale or jump on it.
- ★ DO NOT exceed the maximum weight capacity of 400 lb (180 kg / 28 st).
- ★ DO NOT drop the scale or subject it to strong impacts, as this may damage the sensors.
- ★ DO NOT immerse the scale in water or use chemical cleaning agents. Clean the scale with a slightly damp cloth.
- ★ Pregnant women and children under 14 years old should only use the weight function. Other functions are not intended for them.
- ★ Measurement results from this device are for reference only and not intended to diagnose or treat any disease or abnormalities. Please consult your physician if you have any questions or concerns related to your health.
- ★ Always place the scale on a hard, dry, and flat surface before use.
- ★ Make sure your feet are dry before stepping on the scale.

1. SPECIFICATIONS

Product Name: Smart Body Fat Scale

Model: FI2419

Weight Capacity: 400lb / 180kg / 28st

Increments: 0.1lb / 0.05kg

Weight Units: lb / kg / st

VA Display: 107 × 60mm

Dimensions: 310 × 310 × 35 mm (Handlebar excluded)

Battery: 4 × 1.5V AAA Batteries (inserted)

Platform: Tempered Glass


Auto-off: 15 Seconds

Operation Temperature: 10°C - 40°C

Operation Humidity: 20%- 90%

2. Product Structure

Thank you for choosing this smart body fat scale.


Note:

● How To Wake Up The Scale:

Place it on a flat, hard surface without carpet, and lightly step on the scale to activate it.

(Minimum weight sensing: 5kg)

3. Display Screen

Note: The screen among different models varies, please refer to the actual product for details.

Body Mass Index (BMI) Table

BMI (Suitable for men and women)					
Age	Underweight	Healthy weight	Over weight	Obesity	Severeobesity
18-100	<18.5	18.5~23.9	24~27.9	28~34.9	≥35

Body Mass Index (BMI) = Weight (kg) ÷ Height² (m)


For example: BMI = 70kg ÷ (1.75m* 1.75m) = 22.86

In the range of healthy weight


Note: The BMI index is a commonly used standard to measure body weight and health status.

4. Starfit App

Starfit App can not only track your weight changes and record up to 28 body composition metrics, but also make an analysis and give a health report with the smart body fat scale.

"Measure" Page

More than 20 types of body data

Note: Due to the quick update of the Starfit App, please refer to the actual product for the functions listed above. If there is any change, please update the App to the latest version.

4.1 Install the Starfit App

Search for " Starfit " in App Store or Google play, or scan the QR code below to download the App.

If you use WeChat to scan the QR code, click on the upper right corner on the screen and choose to open it in the browser. App downloaded from other platform or source may be unusable.

Support devices:

- iOS 8.0 and above.
- Android 6.0 and above and equipped with Bluetooth 4.0.
- ※ iPad and tablets are not supported yet.


Scan the QR code
Install the Starfit App

4.2 Register the Starfit app


Use your email address to register a Starfit account or simply log in with your Apple ID or Google account. Fill in your personal information to complete.

Note: Entering incorrect gender, age, or height information may result in inaccurate measurement results


4.3 Connect the Scale with Starfit app


- ① Place the scale on a hard, flat, non-carpeted floor.

- ② Turn on Bluetooth/Network/Location on your phone and remove the plastic insulation strip from the battery compartment to activate the scale.

- ③ Open Starfit, click 'Account' → 'MyDevice' → '+' → 'Searching Bluetooth for pairing' to find the scale 'Runstar-R15' to pair. And now your phone is connected to the scale.


Once connected, the indicator "Ⓢ" will appear on the scale screen, and "Connected" will be shown at the top of the "Measure" page in the app. If "Disconnected" is displayed, it indicates a failed pairing; please repeat the steps above.

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5.1 User Settings

1. Tap SET to activate the scale and keep it on the "Hi" screen.
2. Tap SET again to enter user settings . " USER " will be displayed .

3. Tap<or> to select user code (P0-U1-U2-U3-U4). Then tap SET to set up user details.

* PO is the guest mode. Data is not saved in PO.

4. Setup user details:


- a. Tap<or> to set up age and tap SET to confirm.
- b. Tap<or>to set up gender and tap SET to confirm.
- c. Tap<or>to set up height and tap SET to confirm and save all the user information.

5. Modify user information: Follow the step 1-3 to choose the user you want to modify.

Then long tap SET and follow step 4 to modify the user information.


5.2 Unit Settings

1. Tap SET to activate the scale and keep it on the "Hi" screen.
2. Tap < to select unit (kg, lb, st:lb), Then step onto the scale to weigh.

5.3 Volume Settings

1. Tap SET to activate the scale and keep it on the "Hi" screen.
 2. Tap > to adjust the volume (low , medium , high or mute) .
- Then step onto the scale to weigh.


5.4 Offline Data Review

1. Tap SET to activate the scale and keep it on the "Hi" screen.
2. Long tap<for 3 seconds to enter data review mode.
3. Tap<or> to choose the user you want to review, then tap SET to confirm.
4. Tap<or> to review from DATA 1 to 8.

Each offline user can have a maximum data storage of 8 units.

5. Finally, long tap SET for 3 seconds to exit data review mode.


(Data synced to the app must be reviewed within the app.)

5.5 Delete All Users

1. Tap SET to activate the scale and keep it on the "Hi" screen.
2. Long tap SET for 5 seconds, "DEL" will be displayed. Then tap SET to enter user delete mode.
3. Tap<or> to select "YES" or "NO", then tap SET to confirm. ("YES" to **delete all**, "NO" to cancel.)

*App-synced users can only be deleted via the app.


6. Start Measurement

Before measurement

- **First use:** Remove the plastic insulation strip from the battery compartment to activate the scale.
The screen will display “0.0” and “Hi”.
- **Subsequent use:** Place the scale on a flat, hard surface without carpet. Step on the scale and then step off to activate it. The screen will display “0.0” and “Hi”.

6.1 Correct stance & Correct grip

1. Step onto the scale with bare feet.

Ensure that the soles of your feet are in full contact with the ITO electrodes on the scale (pay attention to left and right placement), as shown in the figure.

2. Pull out the measuring handle. Make sure your thumbs and other fingers are in full contact with the metal electrodes on the handle.

3. Stand upright, naturally extend your arms and keep them straight at a 30°- 45° angle to your body, as shown in the figure.

4. Keep your body stable and avoid touching your abdomen with your arms during measurement.

5. Do not step off the scale until the scale screen or app indicates that the measurement is complete to ensure data accuracy.

Note:

- Measurements must be taken with bare feet, and the skin of both hands and feet must be in full contact with the metal electrode and ITO electrode.
- Place the scale only on a hard, flat surface; do not place any soft materials underneath.
- Remove all packaging from the product to ensure proper skin contact with the electrodes.
- Do not step off the scale until the scale screen or App indicates the measurement is complete.

6.2 Measurement Steps

Please follow the steps below to start the measurement:


1. Step onto the scale and hold the handle firmly with both hands as instructed.
2. After your weight is stabilized, the value will flash three times and lock.
3. Body fat measurement will begin when the colorful dot matrix appears to “dance.”

Please maintain the correct gesture until the measurement is complete.

4. If the scale is not muted, the voice broadcast will play during measurement.

Please follow instructions to complete the process.

5. The scale will automatically shut off after displaying the results for 15 seconds.


The device is turned on.

The weight value flashes three times and locks.


Body fat measurement is begin.

The measurement is complete.


7. Battery and Fault Reminders


Low Battery
(Please replace the batteries.)

Overload (Exceeds maximum
400lb capacity.)

Unstable Weighing
(Please measure again.)

Failed to Obtain Body Fat Data
(Please refer to 8. Frequently Asked
Questions)

8. Frequently Asked Questions

1. The scale failed to connect.

- * Make sure the scale is connected and keep the app on “Measure” page.
- * Ensure that Bluetooth and network permissions for the Starfit APP are enabled. If not, please enable them in the phone settings.
- * Check your smartphone’s operating system version; it must be iOS 8.0 or higher/Android 6.0 or higher.
- * Try deleting the device in the app and reconnecting.

2. No body fat data after weighing.

- * Make sure the scale is connected to the APP via Bluetooth.
Check if the scale is connected to a different phone.
- * Make sure your measuring posture is correct and remain still on the scale for at least 20 seconds until the measurement is complete.
- * Ensure you have entered your height and gender in the Starfit APP.
- * Step onto the scale with bare feet. Make sure your feet and the scale are both dry. If body fat still doesn't show, try applying body lotion to your feet.

3. Incorrect weight data.

- * To ensure accurate measurements, place the scale on a hard and flat surface. Please do not place the scale on carpet.
- * Before use, inspect each foot mat located on the bottom of the scale and ensure they are free of debris or obstructions.
- * If you have not used the scale for a long time, remove the batteries and reinstall them before use.

4. Recalibrate the scale.

- * Remove and reinstall the batteries.
- * Step onto the scale and then step off, wait until the display shows “0.00”.

5. The scale does not turn on.

- * Replace the batteries with new ones.
- * If the scale still does not turn on, please contact the RunStar Customer Service Team for a replacement.

Important Tips:

1. For first-time use, connect the scale to the Starfit app to access all features and measurement data.
If the scale is not connected to the Starfit app, it will only measure weight and BMI.
2. For accurate measurement results, users should remove their shoes and socks and stand barefoot on the scale. Otherwise, the scale will only measure weight and BMI.
3. Do not switch between connected mode and offline mode during measurement. To prevent data loss, we recommend staying in either connected mode or offline mode for the entire measurement process.
4. When multiple users share the scale, make sure the correct user profile is selected in the Starfit app before weighing. Significant differences from the last recorded measurement for the current user may cause the scale to only provide weight and BMI data.

9. FCC Statement

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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