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Warranty Card

User Info:_____	Product Serial No.: _____
Order No.:_____	Product Model: _____
Product Configuration:_____	Purchase Date:_____
Warranty: One Year Warranty Service	Dealer Name:_____

Quick Start Guide

QR Code Prompt
Scan the appropriate QR code for your region and mobile OS to download.
If redirected to the download landing page, select the correct version as instructed on the page.

Charge the ring with the included charging case. Launch the mobile APP, search and connect the ring to begin use.
Read this guide prior to first use, ...h charging, APP download ...ce pairing step by step.

Instructions for Use
Charge the ring by placing it in the charging case. Before first use, charge the ring until it can power on.
While charging, the indicator light on the charging case flashes red.
Once fully charged, the charging case indicator turns off.
After the ring's indicator light starts flashing, open the APP to search for and pair the device.

Additional Notes
The ring may not be found or paired by the APP while it is charging.
If the device has not been used for a long time, charge it before powering on.
Before first use, read this guide and complete charging, APP download and device pairing following the steps.
Ensure the ring is properly connected during charging. The ring's indicator stays solid green when charging starts successfully.

Connect to APP & Pair the Ring
Connect the ring with UBling APP
1. Scan the QR code to download and open the UBling APP.
2. Go to the **Device** page in the APP.
3. Tap **Add Device**.
4. Tap **Bluetooth Search**, locate UBlingR1 and tap to connect.
5. Complete pairing following the on-screen prompts.

Before Pairing, Please Confirm
Bluetooth on your mobile phone is enabled.
The APP has been granted Bluetooth permission and required device permissions.
The ring is charged and removed from charging.
Keep the ring close to your phone; excessive distance may cause search failure.

If the Device Cannot Be Found
If the ring is not detected, check whether the ring has sufficient power, then re-enter the "Add Device" page to search again.
If connection still fails, turn your phone's Bluetooth off and on again, or restart the APP and retry.

Light Status + Wearing Instructions
Light Status

Wearing Instructions
Wear the ring on a comfortable and properly-fitted finger.
Index, middle or ring finger is recommended.
Position the sensor against the underside of your finger for more stable data readings.
Do not wear it too loosely, as this may cause unstable contact between the sensor and skin.

Keep your hand steady during measurement and avoid vigorous shaking.
Data may fluctuate after exercise, when your finger is wet, or if the ring shifts position.
If readings are abnormal, reposition the ring and retake the measurement.

Heart Rate Measurement – Ring indicator flashes green.

SpO2 Measurement – Ring indicator shows solid red.

Charging – Ring indicator flashes red.

Charging Complete – Ring indicator shows solid green.

Plug in the USB cable to charge the charging case. A solid red light indicates charging; it turns solid green once fully charged

When charging the ring, the charging case indicator breathes in green. Once charging completes, the charging case indicator will turn off automatically.

When the charging case power is low (10%), its indicator flashes red once every second

Compliance Statement + Usage Reminders

Usage Reminders

When the battery level is below 5%, the device cannot take measurements. Please charge it before use.
Do not wear the ring or perform measurements while it is charging.
Do not bend, strike, disassemble or attempt to repair the device by yourself.
Keep the sensor area clean to prevent stains from affecting measurement stability.
If the device will not be used for an extended period, charge it before powering on and connecting.

Disclaimer

Measurement results from this product are for daily reference only and are not intended for medical purposes.
Do not use the measurement results for self-diagnosis or treatment.
Consult a doctor promptly if you feel unwell.
Compliance Notice
This device may emit low-power wireless signals during normal operation. Use the device in accordance with this manual.
Do not disassemble or modify the device on your own.
If the device causes interference with other electronic equipment, adjust the operating distance or stop using it and contact after-sales support.

FCC Warning
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and (2) this device muc-cept any interference received, including interference that may cause unde-sired operation.
Any Changes or modifications not expressly approved by the party responsi -ble for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy, if not installed and used in accordance with theinstructions, may cause there is no guarantee that interference will not occur in a particularinstallation, If this equipment does cause harmful interference to radio or television reception,which can be determined by turning the equipment off and on, the user is encouraged to try tocorrect the interference by one or more of the following measures:Reorient or relocate the receiving antenna.-Increase the separation between the equipment and receiver.-Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.Consult the dealer or an experienced radio/TV technician for help The device has been evaluated to meet general RF exposurre requirement. The device can beused in portable exposure condition without restriction.